

# Pirouet™

## Microcentrifuge Tubes

**Pirouet Centrifuge Tubes** are the trusted choice in laboratories all around the world, for applications ranging from 0.5 mL to 50 mL. It all begins with a proprietary, ultra-pure resin that ensures optimal clarity, purity and strength. Plus, over two decades of research and innovation have been instilled into these tubes, resulting in a product that performs to expectation and provides tremendous value.

### Advantages of Pirouet Microcentrifuge Tubes

#### Clarity.

Ultra-Pure, virgin polypropylene delivers crystal clarity and improves visibility.

#### Graduations.


High-contrast graduations facilitate more accurate volume estimation.

#### Sample Recovery.

Select styles feature NoStick® resin technology to ensure optimum sample recovery.

#### Easy to open and close.


Cap design facilitates a simple yet secure seal.



# Pirouet™ 0.5, 1.5, and 2 mL Microcentrifuge Tubes

- 2 cap styles to suit customer preference:
  - Boil-proof caps are ideal for rigorous high-temp applications and tight storage conditions.
  - Large, boil-resistant caps are easy to open & close, with a large labeling surface.
- Crystal-clear polymer facilitates visualization of samples.
- Rated to 20,000 x G.
- Graduations at every 100 µL for easy volume estimation.
- Available with proprietary NoStick resin to ensure maximum sample recovery.
- Certified free of detectable DNase, RNase, DNA, and PCR inhibitors.  
Sterile tubes are also pyrogen-free.

Tubes shown in actual size.



1110-XX  
0.65 mL



1210-XX  
1.7 mL



1260-00  
1.7 mL



1310-XX  
2 mL



## Graduated Microcentrifuge Tubes

Item #	Max Fill	Description	Colors	Packaging
1110-00	0.65 mL	0.5 mL Boil-Proof Cap	-00, -01, -02, -03, -04, -05, -06, -07, -09, -29	1000/unit, 10 units/case
1210-00	1.7 mL	1.5 mL Boil-Proof Cap		500/unit, 10 units/case
1260-00	1.7 mL	1.5 mL Boil-Resistant Large Cap	-00, -01, -02, -03, -04, -05, -06, -09	500/unit, 10 units/case
1310-00	2.0 mL	2.0 mL Boil-Proof Cap	-00, -01, -02, -03, -04, -05, -06, -07, -09, -29	500/unit, 10 units/case

## Graduated Microcentrifuge Tubes, Sterile

1110-S0	0.65 mL	0.5 mL Boil-Proof Cap, Sterile	-00 Clear	1000/unit, 10 units/case
1210-S0	1.7 mL	1.5 mL Boil-Proof Cap, Sterile		500/unit, 10 units/case
1310-S0	2.0 mL	2.0 mL Boil-Proof Cap, Sterile		500/unit, 10 units/case

## NoStick®, Graduated Microcentrifuge Tubes



1110-10	0.65 mL	0.5 mL, NoStick, Boil-Proof Cap	-00 Clear	500/unit, 10 units/case
1210-10	1.7 mL	1.5 mL, NoStick, Boil-Proof Cap	-00 Clear	250/unit, 10 units/case
1310-10	2.0 mL	2.0 mL, NoStick, Boil-Proof Cap	-00 Clear	250/unit, 10 units/case

## Microcentrifuge Tubes without Caps

1270-00	1.7 mL	1.5 mL, Capless	-00 Clear	500/unit, 10 units/case
1370-00	2.0 mL	2.0 mL, Capless	-00 Clear	500/unit, 10 units/case



1270-00  
1.7 mL



1370-00  
2.0 mL



Made in USA